GILSON® pipetman®

DIAMOND TIPS

PIPETMAN® DIAMOND TIPS

Standard is Gilson

Pipetman® Diamond Tips
Pure Performance Tips
for superior results

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HIGHEST QUALITY TIPS FOR SUPERIOR RESULTS

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MADE OF 100% PURE VIRGIN POLYPROPYLENE

No dyes. No additives. No trace metals. Smooth hydrophobic surface prevents liquid retention and ensures complete dispensing.

A PERFECT SEAL

Two solid sealing rings when used with PIPETMAN pipettes ensure more consistent and reliable results each and every time.

ERGONOMIC TIP FITTING

Best match between pipette and tips reduces the strain to fit and eject tips.

GRADUATED VOLUME MARKERS

Volume graduations make it easy to visually check the volume aspirated.

CONTAMINATION PREVENTION WITH PIPETMAN DIAMOND FILTER TIPS

Gilson's high-density polyethylene filter provides an effective barrier against aerosols.

FINE & FLEXIBLE POINT

The last droplet is easier to dispense.

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STERILIZED TIPS

- Sterilized with SAL 10⁻⁶ guaranteed in accordance with ISO 11137.
- Free of detectable RNases, DNases, DNA, RNA, ATP and pyrogens.
- All packs wrapped in medical device quality plastic that provides a complete shield against contamination.
- Sterilization expiration date on each rack.

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
FILTER TIPS

Also available are
PIPETMAN® Filter Tips
offering total protection
against aerosol
contamination and
**Sterilised PIPETMAN®
Tips**, which are gamma
irradiated and certified
free of DNA, RNA,
RNases, DNases, ATP
proteases and trace
metals.

DIAMOND FILTER TIPS Certificate of Conformity

Certificate of Conformity
PIPETMAN® DIAMOND Filter Tips
Contaminant Free Certified

This package contains a high-quality consumable manufactured in a controlled environment following ISO 14644-1. This lot is certified by validation of the manufacturing process following the requirements defined in ISO 18385.

The products are free of the following contaminants:

- Human DNA
- Bacterial DNA
- RNA
- ATP
- DNase
- RNase
- PCR inhibitors
- Endotoxins
- Protein

This lot of tips has been Gamma irradiated for sterilization. The protocol analysis for validation of the sterilization process to ensure a Sterility Assurance Level (SAL) of 10⁻⁶ has been carried out following ISO 11137-2 for establishing the sterilization dosage.

Quality control and subsequent certification is done by independent, accredited laboratories (NF EN ISO/CEI 17025 - COFRAC).

Product certified ISO 18385
Company certified ISO 9001 and ISO 13485

www.gilson.com/ContactUs

Gilson S.A.S.
19 avenue des Entrepreneurs - 95400 Villiers-le-Bel - France
T +33 (0) 1 34 29 50 00 - F +33 (0) 1 34 29 50 20 - info@gilson.com

Hervé Le Dorze
Quality Director

The certification comprises the following tests:

Bio-molecule detection protocols
For all tests, disposable tips are rinsed with sterile DNase/RNase free water called 'liquid extract'. Negative controls, positive and extract controls are made to validate all steps of the different assays.

Human DNA detection < 2 pg
Preparation of sample with sterile and pyrogen free water. Extraction of DNA with a commercial kit. An aliquot of liquid extract is exposed to PCR reaction reagents containing primers (to amplify the β-actin housekeeping gene). A range of human DNA (2, 5, 10, 20 and 50 pg) is also amplified and used as positive control. Amplifications are obtained after cycles of amplification, analysed by electrophoresis on a 2% agarose gel into TBE 0.5X buffer. No DNA have to be detected in the sample.

Bacterial DNA < 2 pg
Preparation of sample with sterile and pyrogen free water. Extraction of DNA with a commercial kit. An aliquot of liquid extract is exposed to PCR reaction reagents containing primer (to amplify RNA 16S gene). A range of bacterial DNA (2, 5, 10, 20 and 50 pg) is also amplified and used as positive control. Amplifications are obtained after cycles of amplification, analysed by electrophoresis on a 2% agarose gel into TBE 0.5X buffer. No DNA have to be detected in the sample.

RNA detection < 1 pg
Preparation of sample with sterile and pyrogen free water. Extraction, purification and RNA concentration are performed with a kit, and RNA converted to cDNA with an other kit, procedure recommended by the manufacturer. The cDNA is exposed to PCR reaction reagents containing primers (to amplify the β-actin housekeeping gene). A range of human RNA (1, 2, 4, 8 and 16 pg) is also amplified and used as positive control. Amplifications are obtained after cycles of amplification, analysed by electrophoresis on a 2% agarose gel into TBE 0.5X buffer. No RNA have to be detected in the sample.

Endotoxin Detection < 0.005 EU/mL
Gilson sterilised DIAMOND Tips certified non-pyrogenic have been tested for bacterial endotoxins. Samples selected at random were tested and validated using the LAL kinetic chromogenic method D with < 0.005 EU/mL sensitivity (European Pharmacopoeia section 2.6.14 methodology for bacterial endotoxin testing).

DNase detection < 12.5 pg
Preparation of sample with sterile and pyrogen free water. 10 µL of liquid extract are incubated 2 hours at 37°C with variable concentrations of RNase plus DNA ladder (0, 10, 50 and 100 ng). Positive controls (0, 10, 50 and 100 ng of DNA + DNase I (12.5 pg) and negative controls (0, 10, 50 and 100 ng of DNA) are incubated in the same conditions. Samples are analysed by electrophoresis (2% agarose gel into TBE 0.5X buffer). The intensity of samples signals is compared with negative and positive controls. The degradation of DNA indicates the presence of DNase in the liquid extract. No degradation have to be observed, the intensity of negative control and the sample have to be identical.

RNase detection < 0.25 ng
Preparation of sample with sterile and pyrogen free water. 10 µL of liquid extract are incubated 10' at room temperature (20-22°C) with 350 ng of RNA (β-2-microglobulin). Positive control (350 ng of RNA + 12.5 ng of RNase A) and negative control (350 ng of RNA) are incubated in the same conditions. Samples are analysed by electrophoresis (2% agarose gel into TBE 0.5X buffer). The intensity of samples signals is compared with negative and positive controls. The degradation of RNA indicates the presence of RNase in the liquid extract. No degradation have to be observed, the intensity of negative control and the sample have to be identical.

ATP detection < 10⁻⁶ mol/L
Preparation of sample with sterile and ATP free water. Evaluation of ATP by luminescence (enzymatic reaction of luciferase (commercial kit)). The intensity of the samples signals is compared with a negative control (DNase/RNase free water) and with ATP standard curve. The intensity of samples signals is compared with negative and positive controls curve. No ATP have to be detected.



Protein detection < 2 µg/mL
Preparation of sample with sterile and pyrogen free water. 100 µL of liquid extract are tested for the detection of proteins by absorbance measurements (commercial kit). Positive control is a standard range supply by the kit, the negative control is sterile and pyrogen free water. Analysis are performed by absorbance reading at the plate reader. The intensity of samples signals is compared with negative and positive controls curve. No protein have to be detected.

Inhibitors of PCR detection
Preparation of sample with sterile and pyrogen free water. Add of 50 pg of control genomic DNA for each sample. Amplification by PCR. Positive control is 50 ng of genomic DNA and an inhibitor standard (e.g. ethanol 0.2%, 0.5%, 1%, 2%, 5%, negative control is sterile and pyrogen free water. A DNA fragment with the same size of amplicon have to be detected.

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

REFILL VERSIONS FOR PIPETMAX® CERTIFIED 384-WELL

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